


BE-GRADIENT BARRIER-FREE



SPECIFICATIONS


TECHNICAL FEATURES

FILLING AND HANDLING

CELL CULTURE

CULTURE MODELS

BE-GRADIENT BARRIER-FREE



BE-Gradient Barrier-Free is our device designed for **electrochemical gradient application to 3D cell cultures** and for creating **barrier models**. It is compatible with any type of **optical microscopy** (inverted phase contrast, confocal, fluorescence, etc). Be-Gradient consists of a cell culture chamber and two channels in direct contact with it. Changing the concentration of one element between the channels on both sides of the central chamber results in the creation of a **chemical gradient** of that element in the central chamber. For barrier and vasculature models, the lateral channels are meant to simulate blood vessels, where endothelial cells can be cultured.

TECHNICAL FEATURES

Be-Gradient Barrier-Free has a horizontal structure for **optimal imaging**, having all structures at the same height. The chamber can confine a hydrogel and so is designed for **3D cell culture**. The two channels flanking the chamber function as **media perfusion** and are directly communicated with the chamber, allowing nutrient and oxygen supply to the central chamber.


For further instructions regarding proper use of the device check the following sections.

CONTENT

The box contains **10 BE-Gradient Barrier-Free** chips that are **individually packaged**. Each chip undergoes sterilization before shipping. It should be stored in a **dry place** which is **not exposed to direct sunlight** at room temperature (**15-25°C**).

	Height	Width	Length	Total volume
Central chamber	200 µm	3 mm	4,5 mm	6,2 µL
Lateral Channels	200 µm	1 mm	50 mm	10 µL
Inlet/Outlet	7 mm	UNF ¼" - 28		130 µL
Medium Reservoir	5 mm	3,6 mm	8,8 mm	185 µL

** The volumes presented on the table are theoretical values calculated for the standard products. Changes in the device features of custom chips may modify the exact channel volume.*




BE-GRADIENT BARRIER-FREE

FILLING AND HANDLING

Before seeding, **prewarm the device** in the incubator overnight to avoid the appearance of air bubbles.

The **chamber must be completely clean and dry** in order to confine the hydrogel in it.

When filling the culture chamber up, place **the tip of the pipette completely vertical into the inlet pinhole**. Inject the liquid with a **continuous and constant flow** to avoid backward flow.




CELL SEEDING/CELL CULTURE

Cell culture without tubing connection


1. Trypsinize and count cells as usual.

2. **Inject the mixture of cells and hydrogel** in the culture chamber by **pipetting through the central inlet** by one of the two following ways:


a. **Inject the solution (10 μ L) slowly** into the chip (A-C). First, the meniscus will enter into the chamber (C). Continue until the chamber is half filled (D,E). Then, inject the remaining solution slowly into the chip through the central outlet until the chamber is fully filled (F-H). When changing the tip from the inlet to the outlet to finish injecting the solution, keep pressing the pipette plunger so air does not enter the pipette tip.




b. **Inject the solution (10 μ L) slowly** into the chip (A). First, the meniscus will enter into the chamber (B,C). Continue slowly until the chamber is fully filled (D-F).



3. **Flip the chip upwards and downwards** every 30 seconds for a few minutes for **homogeneous 3D distribution** during the polymerization of the hydrogel. Disclaimer: the time between flips and incubation might vary according to the type of hydrogel. Check the hydrogel manufacturer's conditions.



4. **Seal the central chamber inlet and outlet** with the adhesive provided. Make sure that the surface is clean to assure proper adhesion.




Important: Note that the channel should never dry up. Take incubation times into account and use the evaporation reservoirs as long as possible.

BE-GRADIENT BARRIER-FREE

CELL CULTURE

Cell culture without tubing connection

5. Add medium to the reservoirs to fill the lateral channels. Gently aspirate the culture medium through the opposite inlet.



a. Add PBS/water to the evaporation reservoir.

b. Additionally, place an **opened Eppendorf's with sterilized water** inside the petri dish to increase humidity and avoid fast evaporation.

6. Cover and incubate until needed.




7. Refresh the medium according to cellular requirements. You may refresh the medium manually at different time points or constantly with a rocker.

Remember that our devices allow the thorough control of the culture conditions and the recreation of hypoxia conditions because COP has very low gas permeability. A correct and frequent replenishing of the media through the lateral channels is important to bring oxygen and nutrients to the culture.

Important: Do not remove the medium inside the channel during medium refreshing. This will prevent the entrance of air bubbles. Thus, when removing the solution, **place the tip of the pipette away from the pinhole** to avoid emptying the channel.

Refreshing medium is required to avoid oxygen and nutrient depletion in the media close to the culture chamber.



Important: Note that the channel should never dry up. Take incubation times into account and use the evaporation reservoirs as long as possible.

FLOW SET UP

Cell culture with tubing connection

Before setting the flow up:

- Before seeding, **sterilize and pre-warm the device and the tubes** overnight at 37 °C.
- On the day of connection, set the system in a laminar flow cabinet.
- **Hydrogel must be fully polymerized** before mounting the perfusion system.
- The channels and inlet/outlet wells should **never be depleted of culture medium**.
- Inlets and outlets are designed to be for **screw connectors (1/4"-28)**.

Connecting the perfusion system

8. Remove the medium from the reservoirs.

9. **Screw the outlet connector** and ensure that the tube is perfectly fixed. **Remove the displaced medium** from the reservoir.

10. **Prime the tube** that will **connect to the inlet before assembling**. This will prevent air bubbles from entering the device.

11. **Screw the inlet connector** and remove the displaced medium from the reservoirs.

13. Once the system is closed, **switch the flow on**.


Check that there are **no leaks** in the system by leaving the pump **running for a couple of minutes** before placing the device in the incubator. Medium reservoirs will be filled up if the connections are not completely closed.

Check our video: Be-Gradient Barrier-Free: Cell Culture


BE-GRADIENT BARRIER-FREE

CULTURE MODELS


1. Cell/spheroid **invasion and migration** by chemotaxis stimuli
2. **Ischaemia model**. Depending on cell type metabolism, cell density, flow conditions, nutrients & oxygen availability, tumour and ischemia models can be performed
3. **Angiogenesis** and **endothelialization**
4. **Blood-brain barrier** models with neurons and glial cells seeded in the central chamber and endothelial cells in the lateral channels




(1)



(2)



(3)



(4)

BEOnChip S.L.
Campus Río Ebro
Mariano Esquillor s/n
Edificio CEMINEM
50018, Zaragoza, Spain

info@beonchip.com
+34 655 16 16 91
www.beonchip.com